

Binding Control Nanoselector Agarose

NBS0BC



Product Type: Alpaca monoclonal, Nanoselector
Host: Alpaca
Applications: IP, Co-IP, Mass Spectrometry, ChIP
Clonality:

Form: Agarose slurry

Hangzhou HuaAn Biotechnology Co.,Ltd.

Orders: 0086-571-88062880

Technical:0086-571-89986345

Service mail: product@huabio.com

www.huabio.com



华安生物
HUABIO

Applications: WB=Western blot IP=Immunoprecipitation IHC=Immunohistochemistry IF=Immunofluorescence FC=Flow cytometry

Note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".



华安生物
HUABIO

Hangzhou HuaAn Biotechnology Co.,Ltd.

Orders: 0086-571-88062880

Technical:0086-571-89986345

Service mail: product@huabio.com

www.huabio.com

Applications: WB=Western blot IP=Immunoprecipitation IHC=Immunohistochemistry IF=Immunofluorescence FC=Flow cytometry