

# HRP Conjugated Anti-DYKDDDDK Tag (FLAG) Antibody

## 0912-3



**Product Type:** Rabbit polyclonal IgG, primary antibodies  
**Species reactivity:** Species independent  
**Applications:** WB

**Description:** The DYKDDDDK peptide has been used extensively as a general epitope tag in expression vectors. Due to its small size, it is unlikely to affect the tagged protein's biochemical properties. This antibody recognizes the DYKDDDDK peptide fused to either the amino- or carboxy-terminus of targeted proteins.

**Conjugate:** HRP-conjugated

**Immunogen:** Synthetic peptide immune sequence is CDYKDDDDK.

**Positive control:** DYKDDDDK Tag recombinant protein

**Recommended Dilutions:**  
**WB** 1:5,000

**Storage Buffer:** 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.

**Storage Instruction:** Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

**Purity:** Immunogen affinity purified.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

 华安生物  
HUABIO  
www.huabio.cn

Images



**Fig1:** Western blot analysis on DYKDDDDK Tag recombinant protein (80kDa).

**Note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物  
HUABIO  
www.huabio.cn