

Images




Fig1: Western blot analysis of Somatostatin on human Somatostatin recombinant protein using anti- Somatostatin antibody at 1/1,000 dilution.




Fig2: Immunohistochemical analysis of paraffin-embedded human lung cancer tissue using anti- Somatostatin antibody. Counter stained with hematoxylin.




Fig3: Immunohistochemical analysis of paraffin-embedded human pancreas tissue using anti- Somatostatin antibody. Counter stained with hematoxylin.

Note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".

Background References

1. Shiraishi M et al. Serum somatostatin in early-stage Parkinson's disease. *Acta Neurol Scand* 121(4):225-9 (2010).
2. Minsel I et al. Somatostatin actions via somatostatin receptors on the ocular surface are modulated by inflammatory processes. *Endocrinology* 150(5):2254-63 (2009).

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物
HUABIO
www.huabio.cn