

Anti-SAS10 Antibody

ER61255



Product Type:	Rabbit polyclonal IgG, primary antibodies
Species reactivity:	Human, Rat, Mouse
Applications:	WB
Molecular Wt:	52 kDa

Description: Essential for gene silencing: has a role in the structure of silenced chromatin. Plays a role in the developing brain. PTM:Phosphorylated upon DNA damage, probably by ATM or ATR. Belongs to the SAS10 family.

Immunogen: Synthesized peptide derived from part region of human protein

Positive control: Amygdala, Colon, Epithelium, Lung, Spleen,

Database links: SwissProt: Q9NQZ2 Human

Recommended Dilutions:
WB 1:500-1:2000

Storage Buffer: Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.

Storage Instruction: Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Purity: Immunogen affinity purified.

Hangzhou HuaAn Biotechnology Co.,Ltd.

Orders: 0086-571-88062880

Technical:0086-571-89986345

Service mail: support@huabio.cn

www.huabio.cn



No Images

Note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".

Hangzhou HuaAn Biotechnology Co.,Ltd.

Orders: 0086-571-88062880

Technical:0086-571-89986345

Service mail: support@huabio.cn

www.huabio.cn



Applications: WB=Western blot IP=Immunoprecipitation IHC=Immunohistochemistry IF=Immunofluorescence FC=Flow cytometry