

Anti-hnRNP CL1 Antibody

ER63694



Product Type:	Rabbit polyclonal IgG, primary antibodies
Species reactivity:	Human
Applications:	WB
Molecular Wt:	32 kDa

Description: May play a role in nucleosome assembly by neutralizing basic proteins such as A and B core hnRNPs. Belongs to the RRM HNRPC family. RALY subfamily. Contains 1 RRM (RNA recognition motif) domain. subcellular location:Component of ribonucleosomes.

Immunogen: The antiserum was produced against synthesized peptide derived from human HNRCL. AA range:31-80

Positive control: Pooled,Uterus,

Database links: SwissProt: O60812 Human

Recommended Dilutions:
WB 1:500-1:2000

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage Instruction: Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Purity: Immunogen affinity purified.

Hangzhou HuaAn Biotechnology Co.,Ltd.

Orders: 0086-571-88062880

Technical:0086-571-89986345

Service mail: support@huabio.cn

www.huabio.cn



Images




Fig1: Western blot analysis of the lysates from RAW264.7 cells using HNRCL antibody.

Note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".