

Anti-Integrin beta 5 Antibody

ER63768



Product Type:	Rabbit polyclonal IgG, primary antibodies
Species reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB
Molecular Wt:	88 kDa

Description: Integrin alpha-V/beta-5 is a receptor for fibronectin. It recognizes the sequence R-G-D in its ligand. Belongs to the integrin beta chain family. Contains 1 PSI domain. Contains 1 VWFA domain. subunit:Heterodimer of an alpha and a beta subunit. Beta-5 associates with alpha-V.

Immunogen: The antiserum was produced against synthesized peptide derived from human Integrin beta5. AA range:721-770

Positive control: Brain,Kidney,Thymic epithelium,

Subcellular location: Membrane.

Database links: SwissProt: P18084 Human

Recommended Dilutions:
WB 1:500-1:2000

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage Instruction: Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Purity: Immunogen affinity purified.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

 华安生物
HUABIO
www.huabio.cn

Images




Fig1: Western Blot analysis of Jurkat cells using Integrin β 5 Polyclonal Antibody diluted at 1:1000

Note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

 华安生物
HUABIO
www.huabio.cn