

# Anti-KOX17 Antibody

ER63816



<b>Product Type:</b>	Rabbit polyclonal IgG, primary antibodies
<b>Species reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC-P
<b>Molecular Wt:</b>	48 kDa

**Description:** May function as a transcription factor. Seems to repress transcription. PTM:Phosphorylated upon DNA damage, probably by ATM or ATR. PTM:Sumoylated. Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 1 SCAN box domain. Contains 4 C2H2-type zinc fingers. tissue specificity:Expressed in many tissues except in heart.

**Immunogen:** The antiserum was produced against synthesized peptide derived from human ZNF24. AA range:10-59

**Positive control:** Bone marrow,Liver,Lung,Lymphoid,Testis,

**Database links:** SwissProt: P17028 Human

**Recommended Dilutions:**

<b>WB</b>	1:500-1:2000
<b>IHC-P</b>	1:100-1:300

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

**Storage Instruction:** Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

**Purity:** Immunogen affinity purified.

Hangzhou HuaAn Biotechnology Co.,Ltd.

Orders: 0086-571-88062880


Technical:0086-571-89986345

Service mail: support@huabio.cn

www.huabio.cn



Images



**Fig1:** Western Blot analysis of HuvEc cells using KOX17 Polyclonal Antibody.

---

**Note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".

---