

Anti-ZIS Antibody

ER65345



| | |
|----------------------------|---|
| Product Type: | Rabbit polyclonal IgG, primary antibodies |
| Species reactivity: | Human, Mouse, Rat |
| Applications: | WB, IHC-P, IF-Cell |
| Molecular Wt: | 34 kDa |

Description: domain:The RanBP2-type zinc fingers mediate binding to RNA. Splice factor required for alternative splicing of TRA2B/SFRS10 transcripts. May interfere with constitutive 5'-splice site selection. PTM:Isoform 2 is phosphorylated on Ser-310 upon DNA damage, probably by ATM or ATR. sequence caution:Intron retention. Belongs to the ZRANB2 family. Contains 2 RanBP2-type zinc fingers. subunit:Interacts with the C-terminal half of SNRNP70, the Arg/Ser-rich domain of XE7 as well as with U2AF1 and CLK1.

Immunogen: The antiserum was produced against synthesized peptide derived from human ZNF265. AA range:106-155

Positive control: Brain,Epithelium,Fetal brain,Fetal kidney,Liver,Skin,Testis,Uterus,

Database links: SwissProt: O95218 Human

Recommended Dilutions:

| | |
|----------------|--------------|
| WB | 1:500-1:2000 |
| IHC-P | 1:100-1:300 |
| IF-cell | 1:100-1:500 |

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage Instruction: Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Purity: Immunogen affinity purified.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

 华安生物
HUABIO
www.huabio.cn

Images




Fig1: Western Blot analysis of Hela cells using ZIS Polyclonal Antibody. Secondary antibody(catalog#:HA1001) was diluted at 1:20,000.

Note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物
HUABIO
www.huabio.cn