

Human LIGHT/TNFSF14 Trimer Protein

HA210223



Product name:	Human LIGHT/TNFSF14 Trimer Protein
Species reactivity:	Human
Protein construction description:	Recombinant Human LIGHT/TNFSF14 Trimer Protein is expressed from HEK293 with His tag and Flag tag at the N-Terminus. It contains Ser89-Val240 Trimer.
Background:	LIGHT, also known as tumor necrosis factor superfamily member 14 (TNFSF14), is a secreted protein of the TNF superfamily. It is recognized by herpesvirus entry mediator (HVEM), as well as decoy receptor 3. This protein has been shown to stimulate the proliferation of T cells, and trigger apoptosis of various tumor cells.
Purity:	> 90% as determined by Tris-Bis PAGE > 95% as determined by HPLC
Endotoxin:	Less than 1EU per µg by the LAL method.
Source:	HEK293
Accession:	O43557-1
Predicted molecular mass:	The protein has a predicted MW of 52.4 kDa. Due to glycosylation, the protein migrates to 53-140 kDa under Non reducing (N) condition based on Tris-Bis PAGE result.
Formulation:	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution:	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage:	Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles. -20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3-6 months after reconstitution. 2-8°C for 2-7 days after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders: 0086-571-88062880

Technical: 0086-571-89986345

Service mail: support@huabio.cn

 华安生物
HUABIO
www.huabio.cn

Applications: WB=Western blot IHC-P=Immunohistochemistry (paraffin) IF-Cell=Immunofluorescence (Cell) IF-Tissue=Immunofluorescence (Tissue) FC=Flow cytometry IP=Immunoprecipitation




Fig1: Human LIGHT Trimer on Tris-Bis PAGE under Non reducing (N) condition. The purity is greater than 90%.




Fig2: The purity of Human LIGHT Trimer is greater than 95% as determined by SEC-HPLC.




Fig3: Immobilized Human LIGHT Trimer, His tag at 1μg/ml (100μl/well) on the plate. Dose response curve for Human HVEM, hFc Tag with the EC50 of 86.3ng/ml determined by ELISA.

Note: All products are “FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE”.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物
HUABIO
www.huabio.cn

Applications:WB=Western blot IHC-P=Immunohistochemistry (paraffin) IF-Cell=Immunofluorescence (Cell) IF-Tissue=Immunofluorescence (Tissue) FC=Flow cytometry IP=Immunoprecipitation