

Rat IL-6, C-His, Flag Tag Protein

HA210629



Product name:	Rat IL-6, C-His, Flag Tag
Species reactivity:	Rat
Bio-Activity:	Testing in progress.
Protein construction description:	A DNA sequence encoding the rat IL-6 protein (P20607) (Phe 25-Thr 211) was expressed with both His, Flag tag at the C-terminus.
Purity:	>95% as determined by SDS-PAGE.
Endotoxin:	Less than 1.0 EU per µg by the LAL method.
Fragment region:	IL-6 (25-211)
Source:	HEK293
Accession:	P20607
Predicted molecular mass:	24.4 kD
Formulation:	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4, 5% Trehalose, 5% mannitol.
Reconstitution:	Reconstitute at 250 µg/ml in sterile water.
Storage:	Please avoid repeated freeze-thaw cycles. Samples are stable for up to twelve months from date of receipt at -20°C to -80°C. It is recommended that aliquot the reconstituted solution to minimize freeze-thaw cycles.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

 华安生物
HUABIO
www.huabio.cn




Fig1: Protein on SDS-PAGE under reducing (R) condition.

Note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".
