

Anti-CD80 Antibody [8-E5-R]

HA601247



Product Type:	Recombinant Mouse monoclonal IgG1, primary antibodies
Species reactivity:	Human
Applications:	WB
Molecular Wt:	Predicted band size: 33 kDa
Clone number:	8-E5-R

Description: CD80/B7-1 is a molecule found on activated B cells and monocytes which provides a costimulatory signal necessary for T cell activation and survival. After engagement of T-cell receptor with antigen in association with major histocompatibility complex class II, a second signal mediated through the binding of B7 to CD28 greatly upregulates the production of multiple lymphokines. B7-1 is a costimulatory molecule for the activation of both CD4+ and CD8+ T lymphocytes that prevents the induction of clonal anergy. Thus, the transfer of B7-1 genes into tumor cells can induce protective immunity and lead to tumor rejection of some tumors in model systems of in vivo tumor growth.

Immunogen: Synthetic peptide within human CD80 aa 74-118 / 288.

Positive control: Human brain tissue lysates.

Subcellular location: Membrane.

Database links: SwissProt: P33681 Human

Recommended Dilutions:

WB 1:1,000

Storage Buffer: PBS (pH7.4), 0.1% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.

Storage Instruction: Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Purity: Protein A affinity purified.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

 华安生物
HUABIO
www.huabio.cn

Images




Fig1: Western blot analysis of CD80 on human brain tissue lysates with Mouse anti-CD80 antibody (HA601247) at 1/1,000 dilution.

Lysates/proteins at 40 µg/Lane.

Predicted band size: 33 kDa

Observed band size: 60 kDa

Exposure time: 24 seconds;

4-20% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (HA601247) at 1/1,000 dilution was used in 5% NFDM/TBST at 4°C overnight. Goat Anti-Mouse IgG - HRP Secondary Antibody (HA1006) at 1/50,000 dilution was used for 1 hour at room temperature.

Note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".

Background References

1. Tekguc M et al. Treg-expressed CTLA-4 depletes CD80/CD86 by trogocytosis, releasing free PD-L1 on antigen-presenting cells. *Proc Natl Acad Sci U S A*. 2021 Jul
2. Sugiura D et al. PD-1 agonism by anti-CD80 inhibits T cell activation and alleviates autoimmunity. *Nat Immunol*. 2022 Mar

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物
HUABIO
www.huabio.cn