

Anti-Junctophilin-2 Antibody [3-2-R] - BSA and Azide free

HA610170



Product Type:	Recombinant Mouse monoclonal IgG1, primary antibodies
Species reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-P
Molecular Wt:	Predicted band size: 74 kDa
Clone number:	3-2-R

Description: Junctophilin-2 (JPH2) is a cardiac specific member of the junctophilins, a newly characterized family of junctional membrane complex proteins important in physically approximating the plasmalemmal L-type calcium channel and the sarcoplasmic reticulum ryanodine receptor for calcium-induced calcium release. It acts probably by anchoring the plasma membrane and endoplasmic reticulum. JPH2 contributes to the construction of skeletal muscle triad junctions, and plays an essential role in heart development.

Immunogen: Synthetic peptide within Human Junctophilin-2 aa 1-50 / 696.

Positive control: Rat skeletal muscle tissue lysates, human skeletal muscle tissue, mouse skeletal muscle tissue, rat skeletal muscle tissue.

Subcellular location: Cell membrane.

Database links: SwissProt: Q9BR39 Human | Q9ET78 Mouse | Q2PS20 Rat

Recommended Dilutions:

WB	1:1,000
IHC-P	1:1,000

Storage Buffer: 1*PBS (pH7.4).

Storage Instruction: Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Purity: Protein A affinity purified.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物
HUABIO
www.huabio.cn

Images




Fig1: Western blot analysis of Junctophilin-2 on rat skeletal muscle tissue lysates with Mouse anti-Junctophilin-2 antibody (HA610170) at 1/1,000 dilution.

Lysates/proteins at 20 µg/Lane.

Predicted band size: 74 kDa

Observed band size: 70 kDa

Exposure time: 3 minutes;

8% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (HA610170) at 1/1,000 dilution was used in 5% NFDM/TBST at 4°C overnight. Goat Anti-Mouse IgG - HRP Secondary Antibody (HA1006) at 1/50,000 dilution was used for 1 hour at room temperature.




Fig2: Immunohistochemical analysis of paraffin-embedded human skeletal muscle tissue with Mouse anti-Junctophilin-2 antibody (HA610170) at 1/1,000 dilution.

The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 20 minutes. The tissues were blocked in 1% BSA for 20 minutes at room temperature, washed with ddH₂O and PBS, and then probed with the primary antibody (HA610170) at 1/1,000 dilution for 1 hour at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.




Fig3: Immunohistochemical analysis of paraffin-embedded mouse skeletal muscle tissue with Mouse anti-Junctophilin-2 antibody (HA610170) at 1/1,000 dilution.

The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 20 minutes. The tissues were blocked in 1% BSA for 20 minutes at room temperature, washed with ddH₂O and PBS, and then probed with the primary antibody (HA610170) at 1/1,000 dilution for 1 hour at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物
HUABIO
www.huabio.cn




Fig4: Immunohistochemical analysis of paraffin-embedded rat skeletal muscle tissue with Mouse anti-Junctophilin-2 antibody (HA610170) at 1/1,000 dilution.

The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 20 minutes. The tissues were blocked in 1% BSA for 20 minutes at room temperature, washed with ddH₂O and PBS, and then probed with the primary antibody (HA610170) at 1/1,000 dilution for 1 hour at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.

Note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".

Background References

1. Takeshima, H., et al. 2000. Junctophilins: a novel family of junctional membrane complex proteins. *Mol. Cell* 6: 11-22.
2. Nishi, M., et al. 2003. Coexpression of Junctophilin type 3 and type 4 in brain. *Brain Res. Mol. Brain Res.* 118: 102-110.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物
HUABIO
www.huabio.cn