

# Anti-ESD Antibody [J0-8-R] - BSA and Azide free

## HA610178



<b>Product Type:</b>	Recombinant Mouse monoclonal IgG1, primary antibodies
<b>Species reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Molecular Wt:</b>	Predicted band size: 31 kDa
<b>Clone number:</b>	J0-8-R

**Description:** ESD (esterase D) is also known as S-formylglutathione hydrolase and is a 282 amino acid protein that is a member of the esterase D family. ESD is highly expressed in placenta, kidney, liver and erythrocytes, and is localized to the cytoplasm, as well as to cytoplasmic vesicles. The main function of ESD is to detoxify formaldehyde while providing energy. Formaldehyde is oxidized by ADH5 which yields S-formylglutathione. ESD then catalyzes the hydrolysis of S-formylglutathione to the reduced forms of formic acid and glutathione. In addition, ESD hydrolyzes a variety of different neutral ester substrates and can act as a carboxylesterase. ESD may also act as a cysteine hydrolase which is inactivated by thiol alkylating agents. ESD gene polymorphism can lead to reduced enzymatic activity which may cause susceptibility to many conditions, including toxic liver cirrhosis, retinoblastoma, obesity and autism.

**Immunogen:** Recombinant protein within Human ESD aa 1-182 / 282.

**Positive control:** K-562 cell lysate, Jurkat cell lysate, HeLa cell lysate, HepG2 cell lysate, SW480 cell lysate, human kidney tissue lysate, mouse kidney tissue lysate, mouse colon tissue lysate, mouse stomach tissue lysate, rat kidney tissue lysate, rat colon tissue lysate, rat stomach tissue lysate.

**Subcellular location:** Cytoplasm. Cytoplasmic vesicle.

**Database links:** SwissProt: P10768 Human | Q9R0P3 Mouse | B0BNE5 Rat

**Recommended Dilutions:**

**WB** 1:1,000-1:20,000

**Storage Buffer:** 1\*PBS (pH7.4).

**Storage Instruction:** Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

**Purity:** Protein A affinity purified.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345


Service mail:support@huabio.cn

华安生物  
HUABIO  
www.huabio.cn

Images

**Fig1:** Western blot analysis of ESD on different lysates with Mouse anti-ESD antibody (HA610178) at 1/1,000 dilution.

- Lane 1: K-562 cell lysate
- Lane 2: Jurkat cell lysate
- Lane 3: HeLa cell lysate
- Lane 4: HepG2 cell lysate
- Lane 5: SW480 cell lysate



Lysates/proteins at 20 µg/Lane.


Predicted band size: 31 kDa  
Observed band size: 31 kDa

Exposure time: 1 minute 13 seconds;  
4-20% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (HA610178) at 1/1,000 dilution was used in 5% NFDM/TBST at 4°C overnight. Goat Anti-Mouse IgG - HRP Secondary Antibody (HA1006) at 1/50,000 dilution was used for 1 hour at room temperature.

**Fig2:** Western blot analysis of ESD on different lysates with Mouse anti-ESD antibody (HA610178) at 1/2,000 dilution.

- Lane 1: Human kidney tissue lysate
- Lane 2: Mouse kidney tissue lysate
- Lane 3: Mouse colon tissue lysate
- Lane 4: Mouse stomach tissue lysate
- Lane 5: Rat kidney tissue lysate
- Lane 6: Rat colon tissue lysate
- Lane 7: Rat stomach tissue lysate



Lysates/proteins at 20 µg/Lane.


Predicted band size: 31 kDa  
Observed band size: 31 kDa

Exposure time: 1 minute 18 seconds;  
4-20% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (HA610178) at 1/2,000 dilution was used in 5% NFDM/TBST at 4°C overnight. Goat Anti-Mouse IgG - HRP Secondary Antibody (HA1006) at 1/50,000 dilution was used for 1 hour at room temperature.

**Fig3:** Western blot analysis of ESD on different lysates with Mouse anti-ESD antibody (HA610178) at 1/20,000 dilution.

Lane 1: HAP1-parental cell lysate  
Lane 2: HAP1-ESD KD cell lysate



Lysates/proteins at 10 µg/Lane.

Predicted band size: 31 kDa

Observed band size: 31 kDa

Exposure time: 10 seconds; ECL: K1801;

4-20% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (HA610178) at 1/20,000 dilution was used in K1803 at 4°C overnight. Goat Anti-Mouse IgG - HRP Secondary Antibody (HA1006) at 1/50,000 dilution was used for 1 hour at room temperature.

**Note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".

## Background References

1. Li W et al. Esterase D enhances type I interferon signal transduction to suppress foot-and-mouth disease virus replication. *Mol Immunol* 75:112-21 (2016).
2. Kaźmierczak M et al. Esterase D and gamma 1 actin level might predict results of induction therapy in patients with acute myeloid leukemia without and with maturation. *Med Oncol* 30(4):725 (2013).

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物  
HUABIO  
www.huabio.cn