

Anti-LXR beta Antibody [PSH0-56]

HA721403



Product Type:	Recombinant Rabbit monoclonal IgG, primary antibodies
Species reactivity:	Human, Mouse, Rat
Applications:	WB, IF-Cell
Molecular Wt:	Predicted band size: 51 kDa
Clone number:	PSH0-56

Description: Liver X receptor beta (LXR- β) is a member of the nuclear receptor family of transcription factors. LXR- β is encoded by the NR1H2 gene (nuclear receptor subfamily 1, group H, member 2). The liver X receptors (LXRs) were originally identified as orphan members of the nuclear receptor superfamily because their ligands were unknown. Like other receptors in the family, LXRs heterodimerize with retinoid X receptor and bind to specific response elements (LXREs) characterized by direct repeats separated by 4 nucleotides. Two genes, alpha (LXRA) and beta, are known to encode LXR proteins. Crystal structure of human liver X receptor β (LXR β) forming heterodimer with its partner retinoid X receptor α (RXR α) on its cognate element, an AGGTCA direct repeat spaced by 4 nt shows an extended X-shaped arrangement, with DNA- and ligand-binding domains crossed. The LXR β core binds DNA via canonical contacts and auxiliary DNA contacts that enhance affinity for the response element. Liver X receptor beta has been shown to interact with NCOA6[8] and Retinoid X receptor alpha.

Immunogen: Synthetic peptide within human LXR beta aa 50-100/460.

Positive control: SW480 cell lysate, HepG2 cell lysate, PC-3M cell lysate, A549 cell lysate, HeLa cell lysate, HT-29 cell lysate, HCT116 cell lysate, HEK-293 cell lysate, mouse brain tissue lysate, mouse liver tissue lysate, rat liver tissue lysate, C2C12, HCT116.

Subcellular location: Nucleus.

Database links: SwissProt: P55055 Human | Q60644 Mouse | Q62755 Rat

Recommended Dilutions:

WB 1:1,000

IF-Cell 1:100

Storage Buffer: 1*TBS (pH7.4), 0.1% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.

Storage Instruction: Shipped at 4°C. Store at +4°C short term (1-2 weeks). It is recommended to aliquot into single-use upon delivery. Store at -20°C long term.

Purity: Protein A affinity purified.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880


Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物
HUABIO
www.huabio.cn

Images

Fig1: Western blot analysis of LXR beta on different lysates with Rabbit anti-LXR beta antibody (HA721403) at 1/1,000 dilution.



Lane 1: SW480 cell lysate (20 µg/Lane)
 Lane 2: HepG2 cell lysate (20 µg/Lane)
 Lane 3: PC-3M cell lysate (20 µg/Lane)
 Lane 4: A549 cell lysate (20 µg/Lane)
 Lane 5: HeLa cell lysate (20 µg/Lane)
 Lane 6: HT-29 cell lysate (20 µg/Lane)
 Lane 7: HCT116 cell lysate (20 µg/Lane)
 Lane 8: HEK-293 cell lysate (20 µg/Lane)
 Lane 9: Mouse brain tissue lysate (40 µg/Lane)
 Lane 10: Mouse liver tissue lysate (40 µg/Lane)
 Lane 11: Rat liver tissue lysate (40 µg/Lane)

Predicted band size: 51 kDa


Observed band size: 51 kDa

Exposure time: 2 minutes;

4-20% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (HA721403) at 1/1,000 dilution was used in 5% NFDM/TBST at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1:100,000 dilution was used for 1 hour at room temperature.

Fig2: Immunocytochemistry analysis of C2C12 cells labeling LXR beta with Rabbit anti-LXR beta antibody (HA721403) at 1/100 dilution.



Cells were fixed in 4% paraformaldehyde for 10 minutes at 37 °C, permeabilized with 0.05% Triton X-100 in PBS for 20 minutes, and then blocked with 2% negative goat serum for 30 minutes at room temperature. Cells were then incubated with Rabbit anti-LXR beta antibody (HA721403) at 1/100 dilution in 2% negative goat serum overnight at 4 °C. Goat Anti-Rabbit IgG H&L (iFluor™ 488, HA1121) was used as the secondary antibody at 1/1,000 dilution. Nuclear DNA was labelled in blue with DAPI.

Beta tubulin (M1305-2, red) was stained at 1/200 dilution overnight at +4 °C. Goat Anti-Mouse IgG H&L (iFluor™ 594, HA1126) was used as the secondary antibody at 1/1,000 dilution.

Hangzhou Huaan Biotechnology Co., Ltd.


Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物
HUABIO
www.huabio.cn

Fig3: Immunocytochemistry analysis of HCT116 cells labeling LXR beta with Rabbit anti-LXR beta antibody (HA721403) at 1/100 dilution.



Cells were fixed in 4% paraformaldehyde for 10 minutes at 37 °C, permeabilized with 0.05% Triton X-100 in PBS for 20 minutes, and then blocked with 2% negative goat serum for 30 minutes at room temperature. Cells were then incubated with Rabbit anti-LXR beta antibody (HA721403) at 1/100 dilution in 2% negative goat serum overnight at 4 °C. Goat Anti-Rabbit IgG H&L (iFluor™ 488, HA1121) was used as the secondary antibody at 1/1,000 dilution. Nuclear DNA was labelled in blue with DAPI.

Beta tubulin (M1305-2, red) was stained at 1/200 dilution overnight at +4 °C. Goat Anti-Mouse IgG H&L (iFluor™ 594, HA1126) was used as the secondary antibody at 1/1,000 dilution.

Fig4: Western blot analysis of LXR beta on different lysates with Rabbit anti-LXR beta antibody (HA721403) at 1/1,000 dilution.

Lane 1: HEK-293-si NT cell lysate

Lane 2: HEK-293-si LXR beta cell lysate


Lysates/proteins at 10 µg/Lane.

Predicted band size: 51 kDa

Observed band size: 51 kDa

Exposure time: 2 minutes;

4-20% SDS-PAGE gel.



Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (HA721403) at 1/1,000 dilution was used in 5% NFDM/TBST at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1:100,000 dilution was used for 1 hour at room temperature.

Note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".

Background References

1. Zhang W et al. The LXRβ-SREBP1 network regulates lipogenic homeostasis by controlling the synthesis of polyunsaturated fatty acids in goat mammary epithelial cells. *J Anim Sci Biotechnol.* 2022 Nov
2. Wolski H et al. Expression of ABCA1 Transporter and LXRA/LXRβ Receptors in Placenta of Women with Late Onset Preeclampsia. *J Clin Med.* 2022 Aug

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

 华安生物
HUAABIO
www.huabio.cn