

Anti-Galectin 9 Antibody [PSH02-45]

HA721822



Product Type:	Recombinant Rabbit monoclonal IgG, primary antibodies
Species reactivity:	Human, Rat
Applications:	WB
Molecular Wt:	Predicted band size: 40 kDa
Clone number:	PSH02-45

Description: The protein has N- and C- terminal carbohydrate-binding domains connected by a link peptide. Multiple alternatively spliced transcript variants have been found for this gene. Galectin-9 is one of the most studied ligands for HAVCR2 (TIM-3) and is expressed on various tumor cells. However, it can also interact with other proteins (CLEC7A, CD137, CD40). For example, an interaction with CD40 on T-cells inhibits their proliferation and induces cell death. Galectin-9 also has important cytoplasmic, intracellular functions and controls AMPK in response to lysosomal damage that can occur upon exposure to endogenous and exogenous membrane damaging agents such as crystalline silica, cholesterol crystals, microbial toxins, proteopathic aggregates such as tau fibrils and amyloids, and signaling pathways inducing lysosomal permeabilization such as those initiated by TRAIL. Mild lysosomal damage, such as that caused by the anti-diabetes drug metformin may contribute to the therapeutic action of metformin by activating AMPK. The mechanism of how Galectin-9 activates AMPK involves recognition of exposed lysosomal luminal glycoproteins such as LAMP1, LAMP2, SCRAMBLIN2, TMEM192, etc., repulsion of deubiquitinating enzyme USP9X, increased K63 ubiquitination of TAK1 (MAP3K7) kinase, which in turn phosphorylates AMPK and activates it. This signaling cascade directly links Galectin-9 intracellular function with ubiquitin systems. Galectin-9, through its regulation of AMPK, a kinase that negatively regulates mTOR, cooperates with Galectin-8-based effects to inactivate mTOR downstream of the lysosomal damaging agents and conditions.

Immunogen: Synthetic peptide within human Galectin 9 aa 51-100 / 355.

Positive control: THP-1 cell lysate, A549 cell lysate, rat liver tissue lysate.

Subcellular location: Cytoplasm, Nucleus, Secreted.

Database links: SwissProt: O00182 Human | P97840 Rat

Recommended Dilutions:

WB 1:2,000

Storage Buffer: PBS (pH7.4), 0.1% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.

Storage Instruction: Shipped at 4°C. Store at +4°C short term (1-2 weeks). It is recommended to aliquot into single-use upon delivery. Store at -20°C long term.

Purity: Protein A affinity purified.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345


Service mail:support@huabio.cn

华安生物
HUABIO
www.huabio.cn

Images

Fig1: Western blot analysis of Galectin 9 on different lysates with Rabbit anti-Galectin 9 antibody (HA721822) at 1/2,000 dilution.

Lane 1: THP-1 cell lysate
Lane 2: A549 cell lysate
Lane 3: Rat liver tissue lysate



Lysates/proteins at 20 µg/Lane.

Predicted band size: 40 kDa
Observed band size: 35 kDa

Exposure time: 2 minutes 18 seconds; ECL: K1802;

4-20% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (HA721822) at 1/2,000 dilution was used in 5% NFDM/TBST at 4°C overnight. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1/50,000 dilution was used for 1 hour at room temperature.

Note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".

Background References

1. Morishita A et al. Galectin-9 in Gastroenterological Cancer. *Int J Mol Sci.* 2023 Mar
2. Yang R et al. Galectin-9 interacts with PD-1 and TIM-3 to regulate T cell death and is a target for cancer immunotherapy. *Nat Commun.* 2021 Feb

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物
HUABIO
www.huabio.cn