

# Anti-WWOX Antibody [PSH15-44]

HA723750



<b>Product Type:</b>	Recombinant Rabbit monoclonal IgG, primary antibodies
<b>Species reactivity:</b>	Human, Mouse, Rat, Monkey
<b>Applications:</b>	WB, IP
<b>Molecular Wt:</b>	Predicted band size: 47 kDa
<b>Clone number:</b>	PSH15-44

**Description:** WW domain-containing oxidoreductase is an enzyme that in humans is encoded by the WWOX gene. WW domain-containing proteins are found in all eukaryotes and play an important role in the regulation of a wide variety of cellular functions such as protein degradation, transcription, and RNA splicing. This gene encodes a protein which contains 2 WW domains and a short-chain dehydrogenase/reductase domain (SRD). The highest normal expression of this gene is detected in hormonally regulated tissues such as testis, ovary, and prostate. This expression pattern and the presence of an SRD domain suggest a role for this gene in steroid metabolism. The encoded protein is more than 90% identical to the mouse protein, which is an essential mediator of tumor necrosis factor-alpha-induced apoptosis, suggesting a similar, important role in apoptosis for the human protein. In addition, there is evidence that this gene behaves as a suppressor of tumor growth. Alternative splicing of this gene generates transcript variants that encode different isoforms. WWOX is also known as human accelerated region 6. It may, therefore, have played a key role in differentiating humans from apes.

**Immunogen:** Recombinant protein within human WWOX aa 1-414.

**Positive control:** MCF7 cell lysate, HepG2 cell lysate, HeLa cell lysate, NIH/3T3 cell lysate, C6 cell lysate, COS-1 cell lysate, Rat brain tissue lysate.

**Subcellular location:** Cytoplasm, Nucleus, Mitochondrion, Golgi apparatus, Lysosome.

**Database links:** SwissProt: Q9NZC7 Human | Q91WL8 Mouse  
Entrez Gene: 292041 Rat

**Recommended Dilutions:**

<b>WB</b>	1:10,000
<b>IP</b>	1-2µg/sample

**Storage Buffer:** 1\*PBS (pH7.4), 0.1% BSA, 40% Glycerol, 0.2% Proclean 950.

**Storage Instruction:** Shipped at 4°C. Store at +4°C short term (1-2 weeks). It is recommended to aliquot into single-use upon delivery. Store at -20°C long term.

**Purity:** Protein A affinity purified.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880


Technical:0086-571-89986345

Service mail:support@huabio.cn

 华安生物  
HUABIO  
www.huabio.cn

## Images

**Fig1:** Western blot analysis of WWOX on different lysates with Rabbit anti-WWOX antibody (HA723750) at 1/10,000 dilution.



Lane 1: MCF7 cell lysate (20 µg/Lane)  
 Lane 2: HepG2 cell lysate (20 µg/Lane)  
 Lane 3: HeLa cell lysate (20 µg/Lane)  
 Lane 4: NIH/3T3 cell lysate (20 µg/Lane)  
 Lane 5: C6 cell lysate (20 µg/Lane)  
 Lane 6: COS-1 cell lysate (20 µg/Lane)  
 Lane 7: Rat brain tissue lysate (30 µg/Lane)


Predicted band size: 47 kDa

Observed band size: 47 kDa

Exposure time: 2 minutes; ECL: K1801;

4-20% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFD/MTBST for 1 hour at room temperature. The primary antibody (HA723750) at 1/10,000 dilution was used in primary antibody dilution (K1803) at 4°C overnight. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1/50,000 dilution was used for 1 hour at room temperature.



**Fig2:** WWOX was immunoprecipitated from 0.2 mg MCF7 cell lysate with HA723750 at 2 µg/10 µl beads. Western blot was performed from the immunoprecipitate using HA723750 at 1/10,000 dilution. HRP Conjugated Anti-Rabbit IgG for IP Nano-secondary antibody at 1/5,000 dilution was used for 1 hour at room temperature.

Lane 1: MCF7 cell lysate (input)  
 Lane 2: HA723750 IP in MCF7 cell lysate  
 Lane 3: Rabbit IgG instead of HA723750 in MCF7 cell lysate

Blocking/Dilution buffer: primary antibody dilution (K1803)

Exposure time: 3 seconds; ECL: K1801

**Note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".

## Background References

- Hsu CY et al. WWOX and Its Binding Proteins in Neurodegeneration. Cells. 2021 Jul
- Baryła I et al. WWOX and metabolic regulation in normal and pathological conditions. J Mol Med (Berl). 2022 Dec

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物  
HUABIO  
www.huabio.cn