

Anti-DXD Antibody [PSH0-73] - BSA and Azide free

HA750652



Product Type:	Recombinant Rabbit monoclonal IgG, primary antibodies
Species reactivity:	Species independent
Applications:	ELISA
Clone number:	PSH0-73

Description: Trastuzumab deruxtecan, sold under the brand name Enhertu, is an antibody-drug conjugate consisting of the humanized monoclonal antibody trastuzumab (Herceptin) covalently linked to the topoisomerase I inhibitor deruxtecan (a derivative of exatecan). It is licensed for the treatment of breast cancer or gastric or gastroesophageal adenocarcinoma. Trastuzumab binds to and blocks signaling through epidermal growth factor receptor 2 (HER2/neu) on cancers that rely on it for growth. Additionally, once bound to HER2 receptors, the antibody is internalized by the cell, carrying the bound deruxtecan along with it, where it interferes with the cell's ability to make DNA structural changes and replicate its DNA during cell division, leading to DNA damage when the cell attempts to replicate itself, destroying the cell. Trastuzumab deruxtecan is the first approved therapy by the US Food and Drug Administration (FDA) targeted to people with the HER2-low breast cancer subtype subset of HER2-negative breast cancer.

Immunogen: DXD-OVA

Recommended Dilutions:
ELISA 1:1,000

Storage Buffer: 1*PBS (pH7.4).

Storage Instruction: Store at +4℃ after thawing. Aliquot store at -20℃. Avoid repeated freeze / thaw cycles.

Purity: Protein A affinity purified.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物
HUABIO
www.huabio.cn

Applications:WB=Western blot IHC-P=Immunohistochemistry (paraffin) IF-Cell=Immunofluorescence (Cell) IF-Tissue=Immunofluorescence (Tissue) FC=Flow cytometry IP=Immunoprecipitation

Images




Fig1: Indirect ELISA analysis of Dxd was performed by coating wells of a 96-well plate with 50 μ L per well of Dxd-SMCC-CCH diluted in carbonate/bicarbonate buffer, at a concentration of 1 μ g/mL overnight at 4 $^{\circ}$ C. Wells of the plate were washed, blocked with 1%BSA blocking buffer, and incubated with 50 μ L per well of Dxd monoclonal antibody serial diluted starting from a concentration of 4 μ g/mL for 1 hour at room temperature. The plate was washed and incubated with 50 μ L per well of an HRP-conjugated goat anti-rabbit IgG secondary antibody at a dilution of 1:35,000 for one hour at room temperature. Detection was performed using an Ultra TMB Substrate for 10 minutes at room temperature in the dark. The reaction was stopped with sulfuric acid and absorbances were read on a spectrophotometer at 450 nm.




Fig2: Competitive ELISA analysis of Dxd was performed by coating wells of a 96-well plate with 50 μ L per well of Dxd-SMCC-CCH diluted in carbonate/bicarbonate buffer, at a concentration of 1 μ g/mL overnight at 4 $^{\circ}$ C. Wells of the plate were washed, blocked with 1%BSA blocking buffer, and incubated with 50 μ L per well of Dxd monoclonal antibody at concentration of 2 μ g/mL with serial diluted Dxd starting from a concentration of 10 μ g/mL for 1 hour at room temperature. The plate was washed and incubated with 50 μ L per well of an HRP-conjugated goat anti-rabbit IgG secondary antibody at a dilution of 1:30,000 for one hour at room temperature. Detection was performed using an Ultra TMB Substrate for 10 minutes at room temperature in the dark. The reaction was stopped with sulfuric acid and absorbances were read on a spectrophotometer at 450 nm.

Note: All products are “FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE”.