



Product Type:	Recombinant Rabbit monoclonal IgG, primary antibodies
Species reactivity:	Human
Applications:	WB, IF-Cell
Molecular Wt:	Predicted band size: 138 kDa
Clone number:	PSH04-03

Description:	HER-2 (also called HER-2/neu, c-erbB2, ERBB2 or neu) is a transmembrane receptor tyrosine kinase. HER-2 is a proto-oncogene, i.e. its activation causes malignant transformation and increases the malignant potential (cell proliferation, invasiveness etc.) of the cells. Amplification of HER-2 gene invariably leads to over-expression of its protein product. The magnitude of over-expression is usually 10-100 folds or even more. Over-expressed HER-2 protein disturbs the HER-receptor family signalling networks, i.e. signalling mediated via EGFR receptor, HER-3 and HER-4. In tumours, HER-2 is over-expressed in 15-25% of primary breast cancers. Metastases usually have the same amplification status as the primary tumours. HER-2 amplification and over-expression are typical features of hormone receptor negative, rapidly growing histologic grade 2-3 tumours. Of the histologic types, Paget's disease is almost invariably HER-2 positive, whereas only a small minority of lobular and tubular carcinomas shows HER-2 amplification. HER-2 amplification and over-expression can also be found in intestinal type gastric and gastroesophageal carcinomas, ovarian carcinomas, high grade endometrial carcinomas and some salivary duct tumours. Low-level copy number increases have been found also in rare cases of lung tumours. Because of its central importance in breast cancer therapy selection, standardization of HER-2 IHC assays and slide interpretation are of outmost clinical and economical importance.
Immunogen:	Synthetic phospho-peptide corresponding to residues surrounding Tyr1248 of Human HER2.
Positive control:	SK-Br-3 cell lysate, SK-Br-3 starved for 4 hours add 200ng/mL EGF for 15 minutes cell lysate, HeLa starved for 4 hours add 200ng/mL EGF for 15 minutes cell lysate, HeLa cells starved for 4 hours add 200ng/mL EGF for 15 minutes.
Subcellular location:	Cell membrane, Cell projection, ruffle membrane; Early endosome, Cytoplasm, perinuclear region, Nucleus.
Database links:	SwissProt: P04626 Human
Recommended Dilutions:	
WB	1:1,000
IF-Cell	1:100
Storage Buffer:	PBS (pH7.4).
Storage Instruction:	Store at +4℃ after thawing. Aliquot store at -20℃. Avoid repeated freeze / thaw cycles.
Purity:	Protein A affinity purified.

Images




Fig1: Western blot analysis of Phospho-HER2 / ErbB2 (Y1248) on different lysates with Rabbit anti-Phospho-HER2 / ErbB2 (Y1248) antibody (HA750926) at 1/1,000 dilution.

Lane 1: SK-Br-3 cell lysate

Lane 2: SK-Br-3 starved for 4 hours add 200ng/mL EGF for 15 minutes cell lysate

Lane 3: HeLa cell lysate

Lane 4: HeLa starved for 4 hours add 200ng/mL EGF for 15 minutes cell lysate

Lane 5: HeLa starved for 4 hours add 200ng/mL EGF for 15 minutes cell lysate, then the membrane treated with λ pp for 1 hour

Lysates/proteins at 20 μ g/Lane.

Predicted band size: 138 kDa


Observed band size: 180 kDa

Exposure time: 24 seconds;

4-20% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (HA750926) at 1/1,000 dilution was used in 5% NFDM/TBST at 4°C overnight. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1/50,000 dilution was used for 1 hour at room temperature.

Fig2: Immunocytochemistry analysis of HeLa cells starved for 4 hours add 200ng/mL EGF for 15 minutes labeling Phospho-HER2 / ErbB2 (Y1248) with Rabbit anti-Phospho-HER2 / ErbB2 (Y1248) antibody (HA750926) at 1/100 dilution.



Cells were fixed in 4% paraformaldehyde for 20 minutes at room temperature, permeabilized with 0.1% Triton X-100 in PBS for 5 minutes at room temperature, then blocked with 1% BSA in 10% negative goat serum for 1 hour at room temperature. Cells were then incubated with Rabbit anti-Phospho-HER2 / ErbB2 (Y1248) antibody (HA750926) at 1/100 dilution in 1% BSA in PBST overnight at 4 °C. Goat Anti-Rabbit IgG H&L (iFluor™ 488, HA1121) was used as the secondary antibody at 1/1,000 dilution. PBS instead of the primary antibody was used as the secondary antibody only control. Nuclear DNA was labelled in blue with DAPI.

Beta tubulin (M1305-2, red) was stained at 1/100 dilution overnight at +4°C. Goat Anti-Mouse IgG H&L (iFluor™ 594, HA1126) was used as the secondary antibody at 1/1,000 dilution.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物
HUABIO
www.huabio.cn

Note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".

Background References

1. Fernandez AI et al. Examination of Low ERBB2 Protein Expression in Breast Cancer Tissue. JAMA Oncol. 2022 Apr
2. Augustin JE et al. Targeting the complexity of ERBB2 biology in gastroesophageal carcinoma. Ann Oncol. 2022 Nov

Hangzhou Huaan Biotechnology Co., Ltd.

Orders: 0086-571-88062880

Technical: 0086-571-89986345

Service mail: support@huabio.cn

 华安生物
HUAABIO
www.huabio.cn

Applications: WB=Western blot IHC-P=Immunohistochemistry (paraffin) IF-Cell=Immunofluorescence (Cell) IF-Tissue=Immunofluorescence (Tissue) FC=Flow cytometry IP=Immunoprecipitation