

Anti-Aquaporin 4 Antibody [PSH06-38] - BSA and Azide free

HA751077



Product Type:	Recombinant Rabbit monoclonal IgG, primary antibodies
Species reactivity:	Human, Mouse, Rat, Cynomolgus monkey, Pig
Applications:	WB, IHC-P, IF-Tissue, IHC-Fr, IP
Molecular Wt:	Predicted band size: 35 kDa
Clone number:	PSH06-38

Description: Aquaporin-4, also known as AQP-4, is a water channel protein encoded by the AQP4 gene in humans. AQP-4 belongs to the aquaporin family of integral membrane proteins that conduct water through the cell membrane. A limited number of aquaporins are found within the central nervous system (CNS): AQP1, 3, 4, 5, 8, 9, and 11, but more exclusive representation of AQP1, 4, and 9 are found in the brain and spinal cord. AQP4 shows the largest presence in the cerebellum and spinal cord grey matter. In the CNS, AQP4 is the most prevalent aquaporin channel, specifically located at the perimicrovessel astrocyte foot processes, glia limitans, and ependyma. In addition, this channel is commonly found facilitating water movement near cerebrospinal fluid and vasculature. Aquaporin-4 was first identified in 1986. It was the first evidence of the existence of water transport channels. The method that was used to discover the existence of the transport channels was through knockout experiments. With this technique they were able to show the significant role of AQP4 in CNS injuries and brain water imbalances. In 1994 the channel was successfully cloned and initially named Mercury-Insensitive Water Channel.

Immunogen: Recombinant protein within human 224-323.

Positive control: Mouse brain tissue lysate, Rat brain tissue lysate, mouse brain tissue, rat brain tissue, human kidney tissue, human brain tissue.

Subcellular location: Cell membrane, Basolateral cell membrane, Endosome membrane, Cell membrane, sarcolemma, Cell projection.

Database links: SwissProt: P55087 Human | P55088 Mouse | P47863 Rat

Recommended Dilutions:

WB	1:2,000
IHC-P	1:2,000-1:10,000
IF-Tissue	1:1,000
IHC-Fr	1:1,000
IP	1-2µg/sample

Storage Buffer: PBS (pH7.4).

Storage Instruction: Store at +4℃ after thawing. Aliquot store at -20℃. Avoid repeated freeze / thaw cycles.

Purity: Protein A affinity purified.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880


Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物
HUABIO
www.huabio.cn

Applications:WB=Western blot IHC-P=Immunohistochemistry (paraffin) IF-Cell=Immunofluorescence (Cell) IF-Tissue=Immunofluorescence (Tissue) FC=Flow cytometry IP=Immunoprecipitation

Images

**Fig1:** Application: IHC-Fr


Species: Mouse

Site: Cerebral cortex

Sample: Frozen section

Antibody concentration: 1:1,000

Antigen retrieval: Not required

**Fig2:** Application: IHC-Fr


Species: Rat

Site: Cerebral cortex

Sample: Frozen section

Antibody concentration: 1:1,000

Antigen retrieval: Not required

**Fig3:** Immunohistochemical analysis of paraffin-embedded human kidney tissue with Rabbit anti-Aquaporin 4 antibody (HA751077) at 1/2,000 dilution and competitor's antibody at 1/2,000 dilution.

The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 20 minutes. The tissues were blocked in 1% BSA for 20 minutes at room temperature, washed with ddH₂O and PBS, and then probed with the primary antibody (HA751077) at 1/2,000 dilution and competitor's antibody at 1/2,000 dilution for 1 hour at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

 华安生物
HUABIO
www.huabio.cn




Fig4: Immunohistochemical analysis of paraffin-embedded human brain tissue with Rabbit anti-Aquaporin 4 antibody (HA751077) at 1/10,000 dilution and competitor's antibody at 1/2,000 dilution.

The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 20 minutes. The tissues were blocked in 1% BSA for 20 minutes at room temperature, washed with ddH₂O and PBS, and then probed with the primary antibody (HA751077) at 1/10,000 dilution and competitor's antibody at 1/2,000 dilution for 1 hour at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.




Fig5: Immunohistochemical analysis of paraffin-embedded mouse brain tissue with Rabbit anti-Aquaporin 4 antibody (HA751077) at 1/10,000 dilution and competitor's antibody at 1/2,000 dilution.

The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 20 minutes. The tissues were blocked in 1% BSA for 20 minutes at room temperature, washed with ddH₂O and PBS, and then probed with the primary antibody (HA751077) at 1/10,000 dilution and competitor's antibody at 1/2,000 dilution for 1 hour at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.




Fig6: Immunohistochemical analysis of paraffin-embedded rat brain tissue with Rabbit anti-Aquaporin 4 antibody (HA751077) at 1/10,000 dilution.

The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 20 minutes. The tissues were blocked in 1% BSA for 20 minutes at room temperature, washed with ddH₂O and PBS, and then probed with the primary antibody (HA751077) at 1/10,000 dilution for 1 hour at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.





Fig7: Immunohistochemical analysis of paraffin-embedded human liver tissue (negative) with Rabbit anti-Aquaporin 4 antibody (HA751077) at 1/10,000 dilution.

The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 20 minutes. The tissues were blocked in 1% BSA for 20 minutes at room temperature, washed with ddH₂O and PBS, and then probed with the primary antibody (HA751077) at 1/10,000 dilution for 1 hour at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.

Fig8: Western blot analysis of Aquaporin 4 on different lysates with Rabbit anti-Aquaporin 4 antibody (HA751077) at 1/2,000 dilution.



Lane 1: Mouse brain tissue lysate (no heat) (40 µg/Lane)
 Lane 2: Mouse liver tissue lysate (negative) (40 µg/Lane)
 Lane 3: Rat brain tissue lysate (no heat) (40 µg/Lane)
 Lane 4: Rat liver tissue lysate (negative) (40 µg/Lane)

Notice: no heat means the lysate is not boiled.

Predicted band size: 35 kDa
 Observed band size: 30 kDa

Exposure time: 16 seconds; ECL: K1801;
 4-20% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (HA751077) at 1/2,000 dilution was used in 5% NFDM/TBST at 4°C overnight. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1/50,000 dilution was used for 1 hour at room temperature.




Fig9: Application: IF-tissue

Species: Mouse

Site: Cerebral cortex

Sample: Paraffin-embedded section

Antibody concentration: 1:1,000

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物
HUAABIO
www.huabio.cn




Fig10: Aquaporin 4 was immunoprecipitated from 0.2 mg mouse brain tissue lysate with HA751077 at 2 μ g/10 μ l beads. Western blot was performed from the immunoprecipitate using HA751077 at 1/1,000 dilution. Mouse anti Rabbit IgG heavy chain (Fc) secondary antibody (M1003-7) at 1/100,000 dilution was used for 1 hour at room temperature.

Lane 1: Mouse brain tissue lysate (input)

Lane 2: HA751077 IP in mouse brain tissue lysate

Lane 3: Rabbit IgG instead of HA751077 in mouse brain tissue lysate

Blocking/Dilution buffer: 5% NFDM/TBST

Exposure time: 1 minute; ECL: K1801

Note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".

Background References

1. Pittock SJ et al. Ravulizumab in Aquaporin-4-Positive Neuromyelitis Optica Spectrum Disorder. *Ann Neurol*. 2023 Jun
2. Kitchen P et al. Targeting Aquaporin-4 Subcellular Localization to Treat Central Nervous System Edema. *Cell*. 2020 May

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物
HUAABIO
www.huabio.cn

Applications:WB=Western blot IHC-P=Immunohistochemistry (paraffin) IF-Cell=Immunofluorescence (Cell) IF-Tissue=Immunofluorescence (Tissue) FC=Flow cytometry IP=Immunoprecipitation