

# HRP Conjugated Anti-6X His tag Antibody

## R1207-3



**Product Type:** Rabbit polyclonal IgG, primary antibodies

**Species reactivity:** Species independent

**Applications:** WB

**Description:** HRP conjugated His-Tag Antibody detects recombinant proteins containing the 6xHis epitope tag. The antibody recognizes the His-tag fused to the amino- or carboxy- termini of targeted proteins in transfected or transformed cells.

**Conjugate:** HRP-conjugated

**Immunogen:** Synthetic peptide immune sequence is HHHHHHHHC.

**Positive control:** C-terminal or N-terminal His-tagged recombinant protein

**Recommended Dilutions:**

**WB** 1:2,000-1:5,000

**Storage Buffer:** 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.

**Storage Instruction:** Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

**Purity:** Immunogen affinity purified.

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880


Technical:0086-571-89986345

Service mail:support@huabio.cn

 华安生物  
HUABIO  
www.huabio.cn

Applications:WB=Western blot IHC-P=Immunohistochemistry (paraffin) IF-Cell=Immunofluorescence (Cell) IF-Tissue=Immunofluorescence (Tissue) FC=Flow cytometry IP=Immunoprecipitation

Images



**Fig1:** Western blot analysis of extracts from cells expressing C-terminal His-tagged protein 6X His-tag antibody:

Lane 1: 1:2,000 diluted antibody

Lane 2: 1:5,000 diluted antibody

**Note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE".

Hangzhou Huaan Biotechnology Co., Ltd.

Orders:0086-571-88062880

Technical:0086-571-89986345

Service mail:support@huabio.cn

华安生物  
HUABIO  
www.huabio.cn